

# Optical Filters: Specifying Filter Spectra

Turan Erdogan, PhD (CTO and Co-founder) Semrock, A Unit of IDEX Corporation


May 31, 2011

www.semrock.com

# **Spectral specifications**

- Semrock always provides

   "manufacturable specifications"
   to customers we guarantee
   the delivered filters will meet
   the specified performance
  - Some filter manufacturers use "design specifications" – these describe what the design will do, and not necessarily what all delivered filters will do
- Example: passband region of exciter and emitter filter pair
  - Note that our knowledge of variations like edge position are taken into account when formulating specifications






# **Spectral specifications**

- Semrock uses a "manufacturable specification" approach to define the bandwidth of a filter passband
- A BrightLine fluorescence filter with a part number FF01-{CWL}/{BW} has a center wavelength of CWL and a guaranteed minimum 93% bandwidth (GMBW) of BW
- Full Width at Half Maximum (FWHM) bandwidth is typically 1% of the CWL wider than the GMBW bandwidth

 $FWHM = GMBW + 0.01 \times CWL$ 

Example – find the FWHM of FF01-520/35 filter (GFP Em):

40 nm = 35 nm + 1% x 520 nm






# **Discrepancy between real and measured data**

- Semrock makes filters with steep edges and deep blocking that "push the limit" of most optical instruments, such as spectrophotometers
- As a result, it can be difficult to measure the filter transmission spectra accurately using established measurement instruments
- There are 3 main discrepancies that appear between the real filter spectrum and the measured spectrum:
  - "Rounding" of sharp spectral features
    - Results from non-zero bandwidth of the spectrophotometer probe beam OR
    - Results from a converging beam used to measure the filter (high cone half-angle)
  - Limited OD measurement range
    - Results from the spectrophotometer "noise floor" detection noise limits how small of a signal can be reliably measured
  - "Sideband measurement artifact" that occurs when measuring filters with a very steep transition from high blocking to high transmission
    - Caused by sidebands of the not perfectly monochromatic probe beam



### **Discrepancy between real and measured data**

Example shows design and measured spectra of an LP02-664RU-25






### How do we handle these data discrepancies?

#### • "Rounding" of sharp spectral features

- We use the finest possible resolution near sharp features and "stitch" the data to more coarse/ efficiently measured data at other wavelengths – never perfect!
- We managed to minimize the cone-half angle of the measurement beam to ~1.0° on Semrock's proprietary spectrophotometer. Note that commercial spectrophotometers have CHAs as large as 2.5° – 3.1°.



### Limited OD measurement range

- Our in-house developed spectrophotometer has improved the range we can measure to (traditionally was OD 5.5). The Near-IR region can be a challenge (requires equipment optimization). Depending on the requirement, it can be costly (in terms of time/ effort) to get precise measurements.
  - Between 320 and 1120 nm, values near or below 3e-7 (OD 6.5) are noise limited
  - For < 320 nm and > 1120 nm, values near or below 3e-6 (OD 5.5) are noise limited

### "Sideband measurement artifact"

 We rely on our knowledge of the spectrophotometer performance and the filter design curve to determine what the actual filter performance is – these measurements are verified periodically using lasers to make singlewavelength OD measurements with a much higher dynamic range (OD values between 8 and 9 can be measured at certain laser wavelengths)





