

Isolation Technologies™ ColumnSaver™

Precolumn Filter

Isolation Technologies ColumnSaver precolumn filters offer all the protection of conventional precolumn filters at 1/3 of the cost. With its convenient direct connect design, changeover time is measured in seconds and no tools or wrenches are needed. The ColumnSaver can also be used as a convenient, fingertight inline filter when used with a 1/16" union.

- ▶ Economical protection for your analytical column
- ▶ Finger tight to 6,000 psi (414 bar) without wrenches or other tools
- ▶ Lower dead volume than conventional filters with holders
- ▶ Universal connection compatible with all column manufacturers' fittings
- ▶ Compatibility with stainless steel or PEEK™ tubing



Ordering Information

ColumnSaver Precolumn

Part Number	Description	Qty.	
9085-05-10	0.5 µm, PEEK	10-pk	\$150.00
9085-20-10	2 µm, PEEK	10-pk	150.00
9085-20T	2 µm, Titanium	ea.	16.00

Eliminate Particulate Contamination

The ColumnSaver PEEK body contains a HiFlow™ filter element and is designed for maximum filtration of particulate matter with minimal dead volume or backpressure. As soon as an increase in backpressure is detected, simply remove and dispose of the ColumnSaver unit and replace it with a new one. The direct-connect design is compatible with all 1/16" 10-32 internal fitting ports from all column manufacturers.

ColumnSaver™, HiFlow™, and Isolation Technologies™ are trademarks of IDEX Health & Science LLC.
PEEK™ is a trademark of Victrex plc.
©2009 IDEX Health & Science LLC

IT1:0-07/2009

IDEX Health & Science LLC
www.idex-hs.com

North/South America +1 866 339 4653 | Europe +49 1801 808 800 | Asia +86 10 6566 9090