

CUSTOMER SUCCESS STORY


Ador Diagnostics Delivered Faster, More Accurate Answers Driven by Multiplexed Amplification and Sensing

The Challenge

- › Multiple reagents: dry and liquid
- › Storage and handling of high and low volume liquid ranges (45 µl to 1,500 µl)
- › Different dry-down conditions for reagents
- › Different assays to be implemented on one platform

Our Solution

- › Separated card manufacturing from reagent functionalization to enable a universal, multi-assay platform
- › Created an on-board storage container with a customized design to meet volume requirements
- › Implemented our Reagent Plug platform for dry reagent storage



“With innovative technologies such as thinXXS’s Reagent Plugs and Wet Reagent Containers, we significantly streamlined reagent handling – a breakthrough that empowered Ador to build a robust and scalable platform.”

– Dr. Eran Zahavy, CEO of Ador Diagnostics