

# Periodic Table of The Analytical Fittings

<b>3750 PSI 1/16</b> <b>M-644-03 / M-650</b> 2-3 in-lbs Headless Micro Knurl PEEK	<b>2500 PSI 1/8</b> <b>P-331 / P-350</b> 3-4 in-lbs Standard Knurl PEEK Tool Required? PEEK	<b>TOOL</b> <b>N-290</b> 3-4 in-lbs Micro Headless Knurl	<b>CONTACT</b> For more information, please visit <a href="http://www.idex-hs.com">www.idex-hs.com</a> or contact: North America CustomerService.hs@idexcorp.com   +1 800 426 0191   +1 360 679 2528 Europe CustomerService.hsEurope@idexcorp.com   +49 1801 808 800   +49 9377 9203-0 Asia CustomerService.hsAsia@idexcorp.com   North Asia +86 10 6566 9090   +86 21 5241 5599 South Asia +91 22 66 43 55 50   Southeast Asia & Oceania +65 6763 6633   Japan +81 48 240 5750	<b>30000 PSI 1/32</b> <b>VHP-700</b> 20 in-lbs 4 mm Hex SS
<b>15000 PSI 1/32</b> <b>UH-904</b> 3-4 in-lbs Headless Knurl PEEK Blend	<b>LEGEND</b> 1 Pressure Rating 2 Tubing Size 3 Fitting / Ferrule Part Number 4 Torque Required 5 Head Style 6 Fitting / Ferrule Material	<b>TOOLS</b> 7 Tool Category 8 Tool Part Number 9 Torque Required 10 Head Style	<b>ICONS</b> Tool Required	<b>NOTE</b> Images are not to scale or shown at actual size.
<b>3750 PSI 1/16</b> <b>M-660 / M-650</b> 2-3 in-lbs Headless Micro Knurl PEEK	<b>4000 PSI .025</b> <b>F-125</b> 3-4 in-lbs Standard Micro Knurl PEEK	<b>5000 PSI 360 µm</b> <b>F-124S</b> 3-4 in-lbs Standard Micro Knurl PEEK	<b>5000 PSI 1/32</b> <b>F-126S</b> 3-4 in-lbs Standard Micro Knurl PEEK	<b>5000 PSI 1/32</b> <b>F-126H</b> 3-4 in-lbs Headless Micro Knurl PEEK
<b>4000 PSI 1/16</b> <b>F-100</b> 3-4 in-lbs Diamond Knurl PCTFE	<b>5000 PSI 1/16</b> <b>F-130</b> 3-4 in-lbs Standard Knurl PEEK	<b>5000 PSI 1/32</b> <b>F-354 / LT-132</b> 4 in-lbs Flush Nut SS / PEEK	<b>5000 PSI 1/32</b> <b>LT-110 / LT-132</b> 4 in-lbs Standard Knurl PEEK	<b>6000 PSI 190 µm</b> <b>F-200 / F-148</b> 3-4 in-lbs Single Wing Delrin / PCTFE
<b>5000 PSI 1/16</b> <b>F-120</b> 3-4 in-lbs Standard Knurl PEEK	<b>5000 PSI 1/16</b> <b>F-127</b> 3-4 in-lbs Standard Knurl PEEK	<b>5000 PSI 1/16</b> <b>F-354 / LT-100</b> 4 in-lbs Flush Nut SS / PEEK	<b>5000 PSI 1/16</b> <b>LT-110 / LT-100</b> 4 in-lbs Standard Knurl PEEK	<b>6000 PSI 190 µm</b> <b>F-300 / F-148</b> 3-4 in-lbs Double Wing PEEK / PCTFE
<b>1000 PSI 1/16</b> <b>P-249</b> 2-3 in-lbs Standard Knurl PEEK	<b>1000 PSI 1/8</b> <b>P-349</b> 2-3 in-lbs Standard Knurl PEEK	<b>1000 PSI 1/8</b> <b>P-331 / P-359</b> 3-4 in-lbs Standard Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> <b>F-356 / P-259</b> 3-4 in-lbs Flush Nut SS / ETFE	<b>1350 PSI 1/16</b> <b>P-255 / P-259</b> 3-4 in-lbs Standard Knurl PEEK / ETFE
<b>1000 PSI 1/8</b> <b>P-100</b> 3-4 in-lbs Diamond Knurl PCTFE	<b>1000 PSI 1/8</b> <b>F-364 / P-359</b> 3-4 in-lbs Flush Nut SS / ETFE	<b>1000 PSI 1/8</b> <b>P-336 / P-359</b> 3-4 in-lbs Headless Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> <b>LT-115 / P-259</b> 3-4 in-lbs Standard Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> <b>P-287 / P-259</b> 3-4 in-lbs Headless Knurl PPS / ETFE
<b>1000 PSI 1/8</b> <b>F-156 / P-359</b> 3-4 in-lbs Female Knurl PEEK / ETFE	<b>1000 PSI 1/8</b> <b>LT-215 / P-359</b> 3-4 in-lbs Standard Knurl PEEK / ETFE	<b>1000 PSI 1/8</b> <b>P-387 / P-359</b> 3-4 in-lbs Headless Knurl PPS / ETFE	<b>1350 PSI 1/16</b> <b>P-232 / P-259</b> 3-4 in-lbs Headless Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> <b>P-420 / P-259</b> 3-4 in-lbs Standard Knurl PEEK / ETFE
<b>1000 PSI 1/16</b> <b>P-229</b> 2-3 in-lbs Standard Knurl PEEK	<b>1000 PSI 1/8</b> <b>P-329</b> 2-3 in-lbs Standard Knurl PEEK	<b>1000 PSI 2.5 mm</b> <b>P-317 / P-366</b> 3-4 in-lbs Standard Knurl PEEK	<b>1350 PSI 1/16</b> <b>P-213 / P-259</b> 3-4 in-lbs Headless Knurl PEEK / ETFE	<b>1350 PSI 1/16</b> <b>P-219 / P-259</b> 3-4 in-lbs Standard Knurl PEEK / ETFE

<b>HAND TIGHTENING/EXTENDER TOOLS</b> <b>F-345</b> 3-4 in-lbs Flush Nut-10-32	<b>F-346</b> 3-4 in-lbs Flush Nut 1/4-20	<b>N-290</b> 3-4 in-lbs Micro Headless Knurl	<b>P-277</b> 3-4 in-lbs Standard Micro Knurl	<b>P-297</b> 3-4 in-lbs Headless Knurl	<b>P-299</b> 3-4 in-lbs Standard Knurl	<b>P-399</b> 3-4 in-lbs Standard Micro Knurl	<b>VHP-1001</b> 10 in-lbs 1/4" Hex	<b>VHP-9000</b> 2-10 in-lbs 6 mm Hex	<b>PRE-SET TORQUE TOOLS</b> <b>VHP-1000</b> 10 in-lbs 1/4" Hex	<b>VHP-2000</b> 16 in-lbs 1/4" Hex	<b>VHP-3000</b> 10 in-lbs 8 mm Hex	<b>P-1000</b> 4 in-lbs Standard Knurl
<b>ADJUSTABLE TORQUE DRIVER AND TORQUE DRIVER EXTENDER TOOLS</b> <b>VHP-4000</b> 2-10 in-lbs 3/16" Drive	<b>N-291</b> 2-10 in-lbs Micro Headless Knurl	<b>P-268</b> 2-10 in-lbs 1/4" Hex	<b>P-278</b> 2-10 in-lbs Female Nut Knurl	<b>P-279</b> 2-10 in-lbs Micro Knurl	<b>P-291</b> 2-10 in-lbs Standard Knurl	<b>P-292</b> 2-10 in-lbs Headless Knurl	<b>F-347</b> 2-10 in-lbs Flush Nut 1/4-28					



Copyright © 2016 IDEX Health & Science LLC. IDEX Health & Science LLC is a Unit of IDEX Corporation. | ID2930

Intelligent Solutions for Life™  
 Fluidics | Optics | Consumables | Assemblies