

# RheBuild® Kit

**5104-999**

RheBuild® Kit for RV500-104 / RV550-104

## Description

The 5104-999 is a maintenance kit for Models RV500-104, and RV550-104 motorized valves. It includes the normal wear parts, tools for disassembly/reassembly and instructions.



## Specifications

### Parts Included in Kit

- (1) Rotor Seal (PEEK)
- (1) Stator Face Seal (PEEK)
- (1) 9/64 Hex Wrench
- Instructions
- **Regulatory REACH and RoHS Compliant: Yes**

**No Dimensional Drawing Required**