

RheBuild® Kit

7116-999

RheBuild® Kit for RV700-116 / RV750-116

Description

The 7116-999 is a maintenance kit for Models RV700-116 and RV750-116 motorized valves. It includes the normal wear parts, tools for disassembly/reassembly and instructions.



Specifications

Parts Included in Kit

- (1) Rotor Seal (PEEK)
- (1) Stator Face Seal (PEEK)
- (1) 9/64 Hex Wrench
- Instructions
- **Regulatory REACH and RoHS Compliant: Yes**

No Dimensional Drawing Required